
GB Model 600 Installation Manual

DE

Model 600 Installationsanleitung

ES

Modelo 600 Manual de instalación

FR

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IT

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JP

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Stannah

GB Introduction

The information contained within this manual has been designed for use by Stannah approved Technicians / Engineers who have the appropriate product training in the following categories:

1. Product Installation
2. Testing and commissioning

The contents of this manual are in the recommended order of placement, which forms the full installation and commissioning of the Stannah 600.

 Pay particular attention to areas marked with this symbol

 To reduce the risk of injury, Personal Protective Equipment should be used where appropriate.

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Contents

Installation Tools.....	4
RR Installation Installation Instructions.....	5
Standard Rail installation Instructions.....	7
Setting options.....	13
Label Positioning.....	20
Battery charger: Positioning Guidance Notes.....	21
Installation Overload Test.....	23
Testing and Commissioning.....	29
Customer Handover.....	47
ECU Controls.....	53
Diagnostics.....	54
Rail end - straight finishing dimensions.....	72
Retractable Rail - Adjustment.....	73
Retractable rail – Bus bar cutting Dimensions.....	74
Electrical Diagrams.....	75

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GB

Installation Tools

- 2.5mm 3mm, 4mm, 5mm, 6mm allen key
- Drill driver
- Tape measure
- Extension bar part # 2992266.
- 13mm Ratchet spanner
- 8mm,10mm Spanner
- Side cutters
- 5 mm, 6 mm Hex key socket
- 13 mm socket
- Torque wrench Range 3Nm—24Nm (2.2lb/ft– 17.6ib/ft)
- PZ3 Bit
- Crimping tool
- Junior Hacksaw
- Spirit level
- Pozi Screwdriver
- Flatblade terminal screwdriver
- Electrical Test meter
- 8mm Twist drill
- 600 rack section

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

Installation Instructions: Retractable Rail

- DE** Montageanleitung: Schiene mit Rückzug.
- ES** Instrucciones de Instalacion: Rail Retráctil.
- FR** Instructions pour l'Installation: Rail Rétractable.
- IT** Istruzioni per l'installazione: Rotaia Retrattile.
- JP** 設置方法: ヒンジレール

- GB**
 1. Assemble first rail leg to suit hand of rail. Note clamp plate.
 - 1a Fit 1x M4 skt button head screw in position 'A'. Fit rail leg to rail in position arrowed and sit rail on stairs.
 2. Raise retractable arm and fit dynamic section. Torque to **12Nm**
 - 2a Lower Retractable rail into closed position to check latching, using 12v battery
 - 2b Fit transit bracket and tighten.
 3. Place 26mm packer under end of rail. Loosen lower rail leg and allow rail to slide down to packer.
 4. Fit 2nd extended rail leg just below rail joint.
 5. Check rail is upright using spirit level to rear of rail and tighten rail leg to rail.
 6. Fit rail joint pins, rail joint plate and fit and secure top rail section.
 7. Fit uppermost rail leg and Lift top of rail to achieve a clearance of **38mm (1½")** between underside of rail and top nosing.

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
- GB**
8. Check rail dimensions:
 - 8a **170mm (6 3/4")** from front top edge of rail to wall.
 - 8b **38mm (1 1/2")** from underside of rail to stair top nosing.
 9. Fix bottom rail leg to stairs with 1 fixing and remove transit bracket.

Ensure appropriate leg support fixing kits are used for the stair tread type, Eg wood fixings for wooden treads.

10. Fit guide arm to upper Retractable rail arm. Note position!
- 10a Fit slider safety assembly to underside of lower rail section and locate guide arm roller into slot.
- 10b Slide complete assembly down to finish level with motor housing and secure.
11. Power the retractable rail up and down using one 12v battery connected to the motor wires to check correct latching. (refer to page 73 if adjustment is required).
12. Fit lower slow/ stop ramp to rear of rail and fit rail end cap.
- 12a Fit lower dynamic cover and refit motor channel cover plate

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


Installation Instructions: Standard Rail

- DE** Montageanleitung: Normschiene
- ES** Instrucciones de Instalacion: Etándar de Ferro
- FR** Instructions pour l'Installation: Rail Standard
- IT** Istruzioni per l'installazione: Ferroviario Normale
- JP** 設置方法：標準レール

- GB**
1. Lay lower section of rail on stairs and slide on rail leg with large halfen to bottom of lower rail section
-  **There are two of these legs. The second leg is to be mounted at the top of the upper rail section. Refer to figure 13a.**
2. Fit lower slowing / stopping ramp and end cap with spring clip. (Position ramp to rear of rail.)
 3. With rail against stair nosings, support rail end 5mm off floor.
 4. Check rail is upright using spirit level on rear of rail.
 5. Tighten rail leg to rail using 13mm ratchet spanner.

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



- GB**
6. Fit and tighten 2nd rail leg to rail below joint. Position rail to achieve:
-  7a **150mm (6")** between wall and front of rail.
 7b **38mm (1½")** between rail and nosing.
 7c **5mm (¼")** between lower end cap and floor.
8. Secure bottom leg with one screw.
 9. Fit joint pins and joint plate to lower section of rail.
 10. Fit top rail section. Refer to page 72, if cutting on site
 11. Fit rail leg with large clamp plate to top of rail for stability.

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

- GB**
12. Push rail joint together. When satisfied with joint, tighten 3 x joint fixings for security and 2 x cone point grub screws for electrical continuity.
 13. Fit further legs as necessary. Ensure:
 - 13a **Max 1.5 Metres (59 1/2") between legs. Max 1.0 Metre (39") between rail joints. Minimum 3 Legs.**
 - 13b **Fit mini leg at top of rail if supplied, for stability**
 14. Screw rail legs to staircase. Ensure rail remains level when securing. If required, fit packers under rail legs to achieve this (part 4303523002 max of 4 per leg).
-  **Ensure appropriate leg support fixings are used for the stair tread type, Eg wood screws for wooden treads**
15. Fit rail leg covers. For Retractable rail, refer to packing list in kit for 1st rail leg cover instructions.

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

GB

16. Remove upper mechanical stop and rack strip stop. Ensure mechanical stop spreader is removed from inside of rail. Remove top section of racking.
 17. Remove carriage cover locking key, carriage covers and upper and lower safety pads.
 - 17a If installing Retractable rail, fit microswitch kit **4301977** to carriage. Refer to fitting instructions within the kit
 18. Using lifting strap, carefully load carriage onto rail.
- Take care not to damage charging brush!**
19. Unclip and lift main ECU. Fit batteries. Switch carriage on via main switch to rear of carriage. Using carriage control buttons, run carriage down rail to a suitable position to fit chair. (Refer to page 53 for control button location.) Switch carriage 'OFF'.

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




- GB**
25. Ensuring terminals are wall facing, connect 24v wire (red) from battery charger.
 - 25a If fitting a retractable rail, feed motor wire to top end of rail under rail legs and secure with spring clips. Cut to length and crimp spade terminals. Connect to terminal block . (Refer to page 76 for connection details)
 26. Slide bus bar back into side extrusion until termination block is positioned against bus bar extrusion. (as photo).
 27. Connect zero volt ring terminal (black wire) from battery charger to underside of rail with shake proof washer.
 28. Cut rack cover to rail length and clip onto rail.
 29. Fit top end cap and secure with spring clip.
 30. If fitting options board, locate and fit 3 spacers and 2 corner fasteners to main carriage ECU tray.
 31. Connect options PCB to main carriage ECU and locate into corner fasteners

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



32

33

- GB**
32. Set handing, CON11
Left Hand installation fit link to pins 1&2
Right Hand installation fit link to pins 2&3.
 33. Set carriage options: SW5
A Call chime
B Travel alarm
C Terminus Alarm
D Retractable rail programming switch. (Refer to page 18)
 34. Set chair options PCB (if fitted): SW7
E Seatbelt switch (Isolation).
F Swivel left allowed
G Swivel right allowed
H Footrest request

Switches E, F, G are for **Starla and Siena chair options** only. For all other chairs, set these switches to the 'OFF' position to ensure correct operation.

Set switch H to the 'ON' position for **Starla / Siena footrest fold button option** and **SL footrest fold button conversion option** only. Set this switch to the 'OFF' position for all other footrest options to ensure correct operation.

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
- GB**
35. Remove Lower interface from packaging and remove interface covers.
 36. Remove carriage lifting strap and fit bolt in top fixing (arrowed), ensuring spacers are located in correct fixing holes.
 37. Fit lower interface onto top bolt.
 - 37a Fit further two fixings (arrowed), level footrest and tighten interface bolts to **24Nm (17.6 lbf/ft)**.
 38. Fit and secure chair to lower interface. Tighten 4 x fixings to **30Nm (22.0 lb/ft)**

Set seat to footrest height to ensure best fit for user.

39. If setting chair to lowest height position, remove break-off tabs from rear cover.

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
- GB**
40. If fitting manual footrest, remove lower nut from bowden cable adjuster, route bowden cable through lower interface and refit nut.
 41. Feed Bowden cable through cable guide, push cable guide onto footrest pivot shaft and locate cable end in footrest plate.
 - 41a Fit and secure clamp.
 42. Adjust Bowden cable
-  **Ensure cable guide is located correctly in footrest top plate (arrowed 42)**
43. Fit footrest safety pad. Hook on at top. Grip as fig 43 and press firmly.
 44. Fit footrest carpet

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

GB

Fitting Upholstery

45. **Siena Chair:** Slide upholstery cover over chair back, observing cable route. Secure with magnets
46. Fit chair base upholstery cover and secure with magnets
47. **Starla Chair:** Connect seat pressure switch wires and fit seat pad.
48. Locate front lip and secure to **8Nm (5.9lb/ft.)**
49. Remove rear casting covers. Fit chair back upholstery and secure to **3Nm (2.2lb/ft)** Replace rear casting covers.

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
GB

Fitting Upholstery

50. **SL Chair** : Fit chair base upholstery and secure fixings to **10Nm (7.4lb/ft)**.
51. Peel double sided tape off rear of infill pad and fit above seat base. (Sofia only)
52. Remove rear cable access covers and fit chair back upholstery. Secure fixings to **10 Nm (7.4lb/ft)**. Replace rear cable access covers

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


GB

53. Remove 2 x link plugs from carriage ECU.
54. Carefully feed chair cables through carriage aperture.
55. Connect to ECU, referring to the following drawings: **Starla /Siena 4309009 page 78 SL 4309019 page 79**
56. Remove charge brush retaining wire.
57. Programme Retractable Rail if fitted as follows
 - With carriage switched 'ON'
 - Ensure switch 'D' on main PCB is in the 'OFF' position. (refer to page 13 for location).
 - Run carriage up and down over entire rail length. NOTE! If powered swivel is fitted, do not allow chair to swivel during this process.
 - Will be displayed on the seven segment display.
 - Continue to operate until display returns to or .
58. If required, rehand keyswitch. **Ensure keyswitch is fitted to upside of chair.**
- 58a Connect key switch and refit all carriage and chair covers.

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


Affix Service Branch contact details :

- DE Befestigen Sie die Kontaktinformationen der Zweigstelle
- ES Fije la información de contacto del servicio técnico
- FR Coller les détails de contact du service après-vente
- IT Applicare i dettagli di contatto della filiale di assistenza
- JP サービス支店の連絡先情報を貼付してください


Carriage Label Position: Left Hand

- DE Position der Etiketten an der Fahreinheit: Linksseitig
- ES Posición de las etiquetas del carro: izquierda
- FR Emplacement des étiquettes de chariot : Côté gauche
- IT Posizione dell'etichetta sul carrello: sinistra
- JP 台座ラベルの位置: 左側


Carriage label Position: Right Hand

- DE Position der Etiketten an der Fahreinheit: Rechtsseitig
- ES Posición de las etiquetas del carro: derecha
- FR Emplacement des étiquettes de chariot : Côté droit
- IT Posizione dell'etichetta sul carrello: destra
- JP 台座ラベルの位置: 右側

Retractable Rail Hazard Label Position R/H Shown

- DE Position der Gefahrenkennzeichnung für die Schiene mit Rückzug auf rechter Seite dargestellt
- ES Posición de la etiqueta de peligro del rail retráctil de la versión de mano derecha
- FR Emplacement de l'étiquette de danger du rail rétractable (côté droit montré)
- IT Posizione dx illustrata dell'etichetta di pericolo della rotaia retrattile
- JP R/H の格納式レール危険標識の位置が示されています

Battery Charger: Position and Connection Guidance Notes

When positioning the battery charger;

- Connection of the battery charger input cable shall be carried out to the latest edition of IEE wiring regulations or the appropriate national standard,
- Fix the battery charger to a suitable surface using fixing screws provided,
- Mechanically protect battery charger input cable for example, inside trunking,
- **Do not** fix battery charger to stairlift rail,
- **Do not** position battery charger in location which obstructs battery charger ventilation slots,
- **Do not** open battery charger enclosure,
- Any work to remove / change the battery charger input plug (if fitted) must be carried out by an appropriately qualified person and meet the relevant national standard,
- **IMPORTANT NOTE : removing / changing the input plug as supplied may invalidate standard approval of battery charger.**
- Avoid positioning the battery charger in location where susceptible to accidental mechanical damage for example, kicked or struck with vacuum cleaner,
- Avoid positioning the battery charger in location where accidental exposure to fluid is likely for example, ground level where floor is cleaned with fluid,
- Avoid positioning the battery charger where it will cause inconvenience to pedestrian users of the staircase.
- When connecting 24v DC cable, cut to correct length. Crimp connections if required as follows: Red +24v = Spade terminal / Black 0v = Ring terminal

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
Overview

It is necessary to complete an important overload safety check upon completion of the installation. A simplified test procedure in accordance with European standard EN81-40 has been created for installed products*.

Details of the test procedure to be done on completion of installation.

1. You should conduct a safe dynamic test by sitting in the chair and riding the stairlift over its full length of travel. Particular attention should be given to rail security and stability.
2. To overload the stairlift a static test should be conducted. You should park the stairlift at the bottom park position, secure your toolbox onto the chair or arm rests, then stand on the front edge of the footrest to create maximum loading, see illustration. Observe the rail security and stability.
3. The same type of static overloading test should be conducted at the upper park position and a mid point on the rail, providing there is a suitable handrail or feature to give you a hand hold while conducting this test at height. Observe the rail security and stability.
4. To carry out a static overload test on a retractable rail, park the stairlift directly over start leg support providing there is a suitable handrail or feature to give you a hand hold. Carry out overloading procedure as above.

What to look for when riding stairlift

The rail legs should remain firmly attached to the floor.
The chair should not tilt excessively.
Expected minimum clearances should be maintained throughout journey.

What to look for during overload test

The rail leg should remain firmly attached to the floor.
No excessive flexing of the rail.

Caution!

- When conducting the overload test ensure the toolbox is secured to the chair to prevent it falling.
- Your toolbox could mark the chair surfaces ensure you adequately protect the chair.
- When conducting the overload test at height ensure you are able to hold something secure to prevent falling.
- When conducting the overload test at height ensure there are no persons directly below you.

*Simplified site test derived from calculations based on rated load and installer weight. (Stannah calculation sheet CS 250-14)

At overload, the product will flex more than usual. The footrest will flex and get closer to the nosings, although it is not expected that it will make contact with the nosings. However, if it does this is acceptable, as in ordinary use the safety pads will ensure that the product stops safely.

A failure is if the overload results in any structural failure of any part of the stairlift or if the fixings pull out, causing a change in the security of any of the rail support fixings. If this occurs, then the stairlift is not suitable for end user use.

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Certificate of Test and Examination after Installation of 600 Model Stairlift to ensure Compliance with CE Directives:

Machinery Directive (2006/42/EC) Low Voltage Directive (2006/95/EC), EMC Directive (2004/108/EC)

Name and Address of User / Installation site: **Tel:**

1. Description

a) **Manufacturer:** *Stannah Stairlifts Ltd*

b) **Stairlift reference number:** **Chair batch No.** **Carriage batch No.**

c) Safe Working Load	24 to 52 degrees = 160kg (352lbs)	Siena / Starla Manual and Powered Swivel
	24 to 52 degrees = 135kg (297lbs)	Sofia / Solus Manual and Powered Swivel

d) Rated speed:	0.1m/s	g) Model type:	600
e) Motor tested on assembly:	YES	h) Type of motor:	DC
f) Control actuator type:	Constant pressure controls	i) Motor protection (type):	Electronic
		j) Type of drive:	Rack and pinion

2. Electrical Tests

a) Electrical supply voltage at time of test:V AC	g) Polarity check:	L-N = L-E = N-E =
b) Battery voltage:V DC	h) 0-volt rail connected and continuity tested:	YES / NO
c) Functional overload protection test:	YES / NO	i) Confirm battery polarity check:	+ ve - ve
d) Tripping time: Seconds	j) AC supply frequency:Hz
e) Control voltage:	24V DC	k) Phase:	Single
f) Mains fuse rating: Amps	l) Motor reversal delay: Seconds

3. Battery Checks

a) Battery charge voltage:	V DC	d) Confirm batteries are fully charged:	YES / NO
b) Confirm battery enclosure is vented:	YES / NO	e) Confirm no leak when battery is fitted:	YES / NO
c) Confirm batteries are securely fitted:	YES / NO	f) Confirm off charge indicator is functional:	YES / NO

4. Mechanical / Operation Test

SWITCH OPERATION			
a) Confirm all on / off stop switches operate correctly:			YES / NO
b) Confirm all control limits operate correctly:			YES / NO
c) Confirm final limits operate correctly:			YES / NO
d) Confirm all sensitive edges / surfaces operate correctly:			YES / NO
e) Confirm underside sensitive surfaces operate correctly and to the appropriate direction of travel			YES / NO
f) Tests to verify the function of the Overspeed Governor and Safety Gear			YES / NO
g) Confirm that all interlock switches operate correctly (e.g on rotating seat, chair arm isolation etc)			YES / NO
h) Confirm retractable rail safety pad switch stops carriage in down direction			YES / NO
i) Confirm that retractable rail safety arm prevents carriage from clashing with the retractable rail.			YES / NO
j) Confirm that retractable rail clutch slips when manually raised			YES / NO
COMPONENT FUNCTIONING			
k) State the stopping distance of the carriage with rated load upon application of the brake:		Up.....mm	Down.....mm
l) Confirm hand winding controls operate correctly.			YES / NO
m) Confirm Over speed Governor has been subjected to factory test for correct tripping speed:			YES / NO
n) Confirm Installation Overload test has been carried out			YES / NO
o) If an alarm is fitted, confirm mode of operation and that the system operates correctly:			YES / NO
p) Confirm clearance of at least 100mm between chair arms and any part of fixed installation whilst in travel position.			YES / NO
q) Ensure minimum 500mm between adjacent stairlifts during full length of travel			YES / NO
r) Confirm that main entrance and exit doors are not obstructed when stairlift not in use			YES / NO
s) Confirm clearance of 1000mm between seat centre and any height restriction during full length of travel			YES / NO
t) Confirm required hazard warning labels are fitted			YES / NO
u) State distance of footrest above floor level: Not to exceed 170mm		Upper..... mm	Lower mm
v) State distance of seat above floor level: Not to exceed 650mm		Upper mm	Lower mm

5. Confirmation

Confirm operating instructions have been handed to the user and the user has been instructed on use of stairlift: YES / NO

6. Declaration

We certify that on (DATE) this stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.

Signed: **Position / Qualifications:**

Introduction

This work method is used for guidance when commissioning a stairlift and completing the Certificate of Test and Examination after the stairlift is installed.

Section 1— Description

- b) Enter details:
 - Contract or order reference number e.g. 99390
 - The chair serial number
 - The carriage serial number
- c) Fit correct weight sticker to carriage label from label pack.
- d) Check that the motor tested label is present on the motor.

Section 2— Electrical Tests

- a) Measure AC supply across the live and neutral terminals at the charger supply input.
- b) Measure battery voltage.
- c) Carry out functional test to confirm current overload device trips effectively.
- d) Electronic trip time (maximum 3 seconds).
- f) Fuse rating (max 2 Amps).
- g) Check polarity at spur or the incoming mains supply.
- h) Carry out a continuity test to ensure the 0 volt wire is mechanically and electrically connected to the rail.
- i) Check polarity of battery connections.
- j) Enter mains AC supply frequency (Hz) for the appropriate country.
- l) Motor reversal delay (min 3 seconds).

Section 3—Battery Tests

- a) Measure battery charge voltage.
- b) Delete as necessary.
- c) Delete as necessary.
- d) Test battery charge condition using the four volt drop principle.
- e) Visually check for any leaks or corrosion.
- f) Check off charge indicator operates.

Section 4 — Mechanical/Operation Test Switch Operation

- a) Delete as necessary.
- b) Delete as necessary.
- c) Delete as necessary.

- d) Delete as necessary.
- e) Delete as necessary.
- f) Test and verify the function of the overspeed governor and safety gear. Refer to page 27 for work method.
- g) All electrical safety devices function correctly.
- h) With downside safety pad removed, press retractable rail safety pad switch whilst placing a downward call. Ensure carriage stops.
- i) With all carriage covers fitted, ensure that the retractable rail safety arm stops the carriage before carriage clashes with the retractable rail. Refer to PDN/ 15/6367
- j) Check that retractable rail clutch allows it to be manually raised.

Component Functioning

- k) Record the stopping distance (must not exceed 20mm).
- l) Check operation of carriage control buttons.
- m) Confirm the tested mark is present on the overspeed governor assembly.
- n) Carry out Installation overload testing. Refer to index on page 3 for work method.
- o) Confirm carriage and Retractable rail alarms and warning LED operate.
- p) Delete as necessary.
- q) Where 2 stairlifts are fitted on the same staircase, ensure that there is a minimum clearance of 500mm between adjacent footrests when unfolded.
- r) Delete as necessary
- s) Delete as necessary
- t) Confirm that hazard warning labels where supplied are fitted to carriage and retractable rail.
- u) Record footrest heights.
- v) Record seat heights

5—Confirmation

The stairlift user must be fully instructed on the use of the stairlift, including the emergency handwinding procedures. All operating instructions should be handed over to the user on completion.

Section 6—Declaration

This section must be completed in fully using the distributor's address and signed.

On completion of the final test and examination and completion of the test certificate, affix the warning labels to the carriage assembly. This is important to ensure compliance with CE directives.

Tests to verify the function of the Overspeed Governor and Safety Gear

Move the carriage to a safe and convenient position on the rail to carry out the work.

Ensure there are no trapping hazards adjacent to the OSG whilst the carriage is moving and the Safety Pad is removed.

- Using a suitable tool, (suggest screwdriver minimum length 100mm or 4”), gently lift the OSG actuating arm whilst driving the stairlift in a downward direction **via the ECU mounted carriage controls buttons** Fig 2, until the OSG switch cuts out the stairlift.
- Check Diagnostics fault code display (9-1 A-6) to prove switch disconnection.
- To reset the OSG, hand-wind the stairlift upwards, using a ¼ inch drive socket on the motor hand winding shaft until the OSG cam plate can be seen to spring back into its travel position Fig 1.
- To reset the ECU, switch the rear mounted carriage switch off and then on.
- Drive the stairlift up about 500mm to completely reset the OSG.
- On completion of the work, test the complete operation of the stairlift.


Figure 1


Figure 2

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Demonstration & Handover

- 1 Ask user to stand clear of the stairlift making sure that they have a clear view of the lift.
- 2 **Getting to know their stairlift** - Talk the user through the key points of the stairlift .
Master power switch / wall switch - Explain the location and use of these switches.
- 3 Inform the customer that the stairlift batteries will not charge with the carriage isolator switch in the '0' position.
- 4 **Power key** - Explain the location and use of the power key. Also explain that they have a spare if the original is lost / damaged..
- 5 **Remote calling controls** - Explain the use of the remote calling controls. Also demonstrate how to replace the batteries.
- 6 **Unfolding** - Demonstrate how to unfold the chair and how to avoid the footrest when lowering the seat.
Using the stairlift - Demonstrate how to get in the seat and put the seatbelt on (if immobiliser seatbelt is fitted explain that the lift will not work until this is on). Also explain the need to keep constant pressure in the direction of travel.
Demonstrate and explain the arm isolation safety cut out, ensure the user is aware that the downside arm has to be fully lowered in order for the stairlift to operate (if fitted).
- 7 Ask the user to sit in the chair, put the seatbelt on and travel up the stairs. Walk up behind the stairlift and explain that it will stop automatically once it reaches the end of the stairs. If an automatic swivel is fitted explain that they need to keep the control pod in the required direction of travel and the chair will swivel.
- 8 **Getting out of the stairlift**- Show the user how to operate the swivel levers (or automatic swivel if fitted) and swivel the chair towards the top of the stairs. Ask them to remove the seatbelt and get out of the chair.
- 9 **Folding** - Demonstrate how to fold the stairlift away.
- 10 **Footrest Fold Button** - (if fitted). Demonstrate footrest fold button option.
- 11 Now ask the user to use the stairlift to go down the stairs, covering any of the above that the user seems unsure about.
- 12 Explain the use of the hand-winding controls (including the 2-way automatic swivel if fitted).
- 13 Explain the location of the safety edges and give a demonstration on how they work.
- 14 Check that all of the above has been understood. If they are unsure on any of the above go through it with them again.
- 15 Complete handover document.
- 16 **2-Way Automatic swivel** - (if fitted) Demonstrate how this is used explaining that it will only work when at the bottom of the stairs.
- 17 **2-Way Manual swivel** - (if fitted) Demonstrate how this is used explaining that it will only work when at the bottom of the stairs.
- 18 **Retractable Rail**. Demonstrate how the retractable rail operates, explaining that the area around the retractable rail should be clear before moving the stairlift.
- 19

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
GB

Main ECU controls

- LED9-** **Status Display** - Operator basic trouble shooting guide
- SW1-** **Carriage control button up/down**
SW2- **Carriage control button up/down**
 These buttons allow the carriage to be operated at half speed.
- SW3-** Latch switch for normal chair operation and to display user faults (0-9) on fault display LED.
 De-latch switch to activate carriage control buttons and to enable engineer fault diagnostic mode.
- SW4-** **Mode Advance button** changes Status display LED (LED9)
- LED's 1- 8 Fault Display.** Illuminate to Identify fault on product.

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

The **Status display** is visible with the carriage cover in place and will show the following basic user fault codes:

- Charge voltage not present.
- 1. Battery voltage too low.
- 2. Ready for operation.
- 3. Safety pad obstructed.
- 4. Chair-arm rocker switch or key switch is OFF.
- 5. Downside chair arm not lowered or seat not loaded.
- 6. Isolator seatbelt not worn. (If Fitted)
- 7. Seat not centred.
- 8. Switch 3 on ECU not latched.
- H-h Switch 3 on ECU not latched. (A stored fault is present, Refer to page 56).
- 9. Fault present. Refer to fault codes.
- F. Software fault. Replace main ECU.


Example Fault:

- 1 Stairlift has stopped. Customer reports '9' is displayed in the **Status Display**. *Status Display showing '9' 'Fault present'*.

- 2 Remove Carriage Cover to access main ECU and **Fault Display** LED's. Press and de-latch SW 3 on main ECU to enable fault mode.


Status Display continues to show '9' on main ECU.

Status Display	FAULT
9	A
Fault Group 3	1 No power supply detected to option ECU.
	2 Battery voltage in excess of 30V.
	3 Current limit exceeded 24A Instantaneous
	4 Current limit exceeded 16A for 1 second
	5 Powered swivel motor stalled (Starla or Saxon) Powered swivel or footrest motor current limit exceeded

Fault Display LED's show fault group '3'. Fault group '3' indicates that it will be one of the faults from 1 to 6.

- 3 For specific fault information, press and hold **Mode Advance** button.

Status display will change from '9' to 'A'.

Status Display	FAULT
9	A
Fault Group 3	1 No power supply detected to option ECU.
	2 Battery voltage in excess of 30V.
	3 Current limit exceeded 24A Instantaneous
	4 Current limit exceeded 16A for 1 second Powered swivel motor stalled (Starla or Saxon) Powered swivel or footrest motor current limit exceeded

Fault display changes, to show specific fault. In this case '4'. "Current limit exceeded 16A for 1 second"

Fault display shows specific fault '4'.

If '-', '1' or '2' is shown in the **Status Display** window, removing the carriage cover and delatching SW3 will show **Fault Display LED's** illuminated to indicate which options are fitted to the stairlift.

Status Display '1' or '2'		
Fault Display LED's	LED 1 On } Saxon chair On } Starla/Siena chair Off } SL chair	LED 5 On - Power swivel left allowed
	LED 2 On } Saxon chair Off } Starla/Siena chair On }	LED 6 On - Power swivel right allowed
	LED 3 On - Powered footrest present	LED 7 On = L/H stairlift Off = R/H stairlift
	LED 4 On - Manual 2 way swivel.	LED 8 On - Option PCB fitted

If '3' is shown in the **Status Display** window, (**Safety pad obstructed**), removing the carriage cover and delatching SW 3 will show **Fault Display LED's** illuminated to indicate which safety pad may be obstructed on the stairlift.

Status Display '3'		
Fault Display LED's	LED 1 On Bottom footrest safety pad activated	LED 5 On Bottom interface safety pad activated
	LED 2 <i>Not Assigned</i>	LED 6 <i>Not Assigned</i>
	LED 3 On Right footrest/interface safety pad activated	LED 7 On Left footrest /interface safety pad activated
	LED 4 On Right carriage safety pad activated	LED 8 On Left carriage safety pad activated.

If '9' is displayed in the '**Status Display**' window, refer to the following chart to identify the fault.

Fault Codes

Fault Display LED's	Status Display		FAULT	ACTION
	9	A		
Fault Group 1	1		Battery voltage below 23V. Carriage will operate from carriage control buttons only.	Check Charger output / Replace Batteries
	2		Battery voltage below 21V (Carriage will NOT operate from any controls)	Check Charger output / Replace Batteries
	5		Ultimate switch activated	Check for over travel / switch / wiring
	6		OSG switch activated	Check OSG / switch / wiring loom
	7		Ultimate Line switch activated	Check chair arm/ OSG / ultimate switches or looms. Identified ULT 1- 4
Fault Group 2	1		Left and Right slowing switches activated simultaneously	Check switches / wiring loom
	2		Stop Switch activated before slowing switch	Check switch operation / connections
	3		Dual calls from chair arm controls	Check arm control / wiring loom
	4		Dual calls from remote controls	Check IR PCB connections
	5		SL Chair: Seat raised and loaded simultaneously	Check seat switches / loom
	6		'A' Circuit shorted to ground. (Chassis)	Identify and check ' A ' wiring looms
	7		'B' Circuit shorted to ground. (Chassis)	Identify and check ' B ' wiring looms
Fault Group 3	1		No power supply detected to option PCB.	Check fuse F2 (3.15A)
	2		Battery voltage in excess of 30V.	Check charger output / motor
	3		Current limit exceeded 24A Instantaneous	Excess load / Check Motor
	4		Current limit exceeded 16A for 1 second	Excess load / Check Motor
	5		Powered swivel motor stalled	Excess load / Check connections
	6		Powered swivel or footrest motor current limit exceeded	Excess load. Check options motors
Fault Group 4	1		Brake not detected	Check brake and brake connections
	2		Main motor not detected.	Check main motor and connections
	3		Footrest motor not detected	Check footrest motor and connections
	4		Swivel motor 1 not detected	Check swivel motor and connections
	5		Swivel motor 2 not detected	Check swivel motor and connections
	6		Main motor current drops below 1A when travelling upwards	Check motor connections / Check handing plug if travelling downwards

Fault Codes

Fault Display LED's

Status Display		FAULT		ACTION	
9	A				
Fault Group 5	1	Starla manual 2 way swivel interlock motor timed out	Check change over mechanism		
	2	Starla manual 2 way swivel. Both switches activated Simultaneously.	Check switches / check wiring		
	3	Starla / Siena power swivel encoder fault. All outputs 'off	Check encoder / connections/ slot alignment		
	4	Starla / Siena power swivel centred switch fault. Switch still activated when encoder centred	Check switch / encoder connections		
	5	Starla power swivel encoder sequence error.	Check / replace encoder		
	6	Starla power swivel encoder incorrect sequence in travel	Refer to instruction below to reprogramme encoder / Replace encoder		
Fault Group 6	1	Main ECU failure	Replace main ECU		
	2	Option PCB failure	Replace option PCB		
	3	Intermittent connection between main ECU and option PCB	Check connection between boards		
Fault Group 7	1	Con 7 wiring loom short circuit between A&B CCTS	Check wiring loom con 7		
	2	Con 9 wiring loom short circuit between A&B CCTS	Check wiring loom con 9		
	3	Con 10 wiring loom short circuit between A&B CCTS	Check wiring loom con 10		
	4	Con 17 wiring loom short circuit between A&B CCTS	Check wiring loom con 17		
	5	Con 18 wiring loom short circuit between A&B CCTS	Check wiring loom con 18		
Fault Group 8	2	Stairlift handing not specified.	Check handing plug fitted correctly		
	3	Starla / Saxon swivel wiring incorrect.	Check chair loom con10		
	4	New SL swivel detected (SL encoder) not supported.	Replace main ECU to the latest version		
F	F	Main ECU software failure	Replace Main ECU		

* If Fault Code identifies switch / wiring fault and testing the circuit proves no fault. ECU failure should be assumed.

Resetting Encoder following a sequence fault

If Resetting the power does not clear an encoder sequence fault, this may be due to the fault occurring during an inwards travel operation. To reset the fault use the following method:

1. Switch power 'OFF'.
2. Set switches 'F' and 'G' on SW7 on options board to the 'OFF' position. (Refer to page 9 for switch location).
3. Switch power 'ON'
4. Switch power 'OFF' and reset switches to desired configuration (Refer to page 9 for switch configuration).
5. Switch power 'ON'.

Accessing Stored Fault Codes

To access this feature:


1. With carriage switched on
2. Remove carriage cover to access the main ECU
3. With status display showing '2', De-latch SW3
4. If there is a stored fault code, the Status Display (LED 9) will cycle between 'H' and 'h'. The Fault Display LED's will cycle simultaneously.

Example: H4-h3 Indicates '*Footrest motor not detected*' as the last stored fault. (Equivalent to fault 9,4-A,3).

5. To clear the stored fault, press and hold SW4 for 4 seconds. The stored fault will be permanently removed.
6. Re-latch SW3 to return stairlift to normal operation.

NOTE. A stored fault cannot be displayed if the Status Display is indicating '9'.

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




Rail end – Straight finishing dimensions

- ES** Final de carril/riel – Dimensiones de una terminacion de carril/riel recta
- IT** Ferroviario fine – Dritto finale dimension
- DE** Ende der Schiene – Finish dimension

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Retractable Rail Adjustment

GB Loosen screw 'A' before adjustment of screw 'B'. After adjustment, re-torque screw 'A' to 12NM.

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






Retractable rail – Bus bar cutting Dimensions

- DE Schiene mit Rückzug – Schneidmaße der Stromschiene
- ES Raíl retráctil – Longitud de corte de las barras de distribución
- FR Rail rétractable – Dimensions de coupe de barre omnibus
- IT Rotaia retrattile – Dimensioni di taglio delle barre collettrici
- JP ヒンジレールバス・バー切断寸法



600 RETRACTABLE RAIL CONNECTION DIAGRAM

IMPORTANT NOTES :


- BATTERY CHARGER MUST NOT BE MOUNTED TO THE RAIL
- MAINS INPUT CABLE MUST NOT BE SECURED TO OR STORED INSIDE THE RAIL

STAIRLIFT SUPPLY WIRING MUST BE CONNECTED IN ACCORDANCE WITH THE APPROPRIATE NATIONAL STANDARD AND THE STANNAH STAIRLIFT INSTALLATION TRAINING PROCESS.


600 BUS BAR BATTERY CHARGER CONNECTION DIAGRAM

RIGHT HAND


LEFT HAND


600 RETRACTABLE RAIL POWER LOOM ROUTING

FIGURE 1


600 BUS BAR BATTERY CHARGER & RETRACTABLE RAIL CONNECTION DIAGRAM


600 DEALER BUS BAR BATTERY CHARGER CONN DIAG

RIGHT HAND

LEFT HAND

600 RETRACTABLE RAIL POWER LOOM ROUTING



FIGURE 1


600 DEALER BUS BAR BATTERY CHARGER & RETRACTABLE RAIL CONN DIAG


FOR DETAILS OF BATTERY CHARGER CONNECTION REFER TO DIAGRAM 4309210.

INTERFACE & FOOTREST

SL CHAIR: KEYSWITCH NOT REQUIRED
CONNECT WIRES USING 2609158


SL CHAIR


The circuits in this system are in two groups, 'A' and 'B'.
The two groups use a 24VDC pulsing signal.
The 'A' pulses are timed differently to the 'B' pulses.
This aids the electronics detect incorrect wiring and faults.

CARRIAGE


SL 2 WAY POWERED SWIVEL



FOR 1 WAY SL POWERED SWIVEL ONLY CONNECT 1 MOTOR (LEFT OR RIGHT SWIVEL)


RETRACTABLE RAIL


POWERED FOOTREST


INTERFACE & FOOTREST


SIENA WIDE ARM CHAIR



2 WAY POWERED SWIVEL

CARRIAGE






POWERED FOOTREST





The circuits in this system are in two groups, 'A' and 'B'. The two groups use a 24VDC pulsing signal. The 'A' pulses are timed differently to the 'B' pulses. This aids the electronics detect incorrect wiring and faults.

600 SIENA WIDE ARM CHAIR CONNECTION DIAGRAM 4309028 C





KEY
 (M) = MAIN CONTROL PCB
 (O) = OPTIONS PCB
 (C) = CHAIR ARM PCB
 (E) = ENCODER TRANSMITTER PCB
 (R) = ENCODER RECEIVER PCB
 (IR) = INFRA RED RECEIVER
 CON = CONNECTOR
 M = MOTOR
 B = BRAKE

KEY

- (M) = MAIN CONTROL PCB
- (O) = OPTIONS PCB
- CON = CONNECTOR
 MOTOR
 BRAKE

600 WITH SL CHAIR KEY WIRING DIAGRAM 4309022C

KEY

- (M) = MAIN CONTROL PCB
- (O) = OPTIONS PCB
- (C) = CHAIR ARM PCB
- (ET) = ENCODER TRANSMITTER PCB
- (IR) = INFRA RED RECEIVER
- CON = CONNECTOR
- ⊕ MOTOR
- ⊖ BRAKE

'Original instructions in English'