GB		Page
	General Installation Instructions	2
	Electrical Wiring Diagrams	11
	Test and inspection methods	21
	Test and inspection methods - Overspeed Governor	27
	Demonstration and handover	28
	Tools and equipment	30
	Warning Labels	32
	Unconfigured Rails	33

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


- GB**
1. Place first section of rail on stairs.
Fit slowing and stopping ramps to bottom of rail.
 2. Fit bottom end cap.
 3. Fit legs to first section of rail.
 4. Ensuring you have the correct distance from floor, nosing and wall, secure bottom leg with one screw.

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
GB

5. Fit joint plate and slide joint plate in.
6. Fit second section of rail. Fit top leg for stability and joint alignment. When satisfied with joint, tighten fixings and fit further legs as necessary. Fix legs to staircase. Ensure rail is level when securing. Ensure 2 metre maximum distance between supports.
7. Remove mechanical stop and top section of rack.
8. Remove carriage from packaging. Remove support strap and fixings. Remove top cover as shown.

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
11

GB

9. Remove carriage safety pads.
 - a) The fixing on both sides of safety pad must be unclipped before lifting safety pad off carriage.
10. Fit trailing cable to carriage by threading connector through slot on underside of carriage, ensuring enough cable is passed through to reach ECU and secure with cable clip.
 - b) Route cable as shown for left hand lift.
 - c) Route cable as shown for right hand lift.
11. Slide carriage onto rail.

Use safety manual handling procedures when lifting carriage.

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
GB

12. Fit batteries and wire in order as shown.
 - 1 - (Red)
 - 2 - (Black - Fuse Link)
 - 3 - (Red - Fuse Link)
 - 4 - (Black)
13. Refit safety pads ensuring correct orientation. Run carriage down to bottom of rail using hand control and link plugs.
14. Adjust bogie block to obtain light clamping force, recommended torque 3Nm (2.2 lbf ft).
15.  Lock bogie block by securing two side fixings, mandatory torque 24Nm (17.7 lbf ft).

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


19

GB

16. Decide on the hole configuration to be used for chair mounting, fit spacers and socket set screw into carriage side.
- a) Socket set screw
 - b) Spacers
17.  Fit chair using socket set screw for alignment. Tighten screws to mandatory torque 24Nm (17.7 lbf ft), ensuring the chair is level. Replace socket set screw with cap head screw and torque to 24Nm (17.7 lbf ft).
- a) Socket set screw
 - b) Feed looms from carriage through interface.
18. Fit rear strap using fixings supplied.
19.  Ensure that spacer is fitted behind strap.

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GB

20. Plug connectors in and refit covers.
21. Fit upholstery using fixings supplied. Also, if seat belt option has been specified, fit using instructions supplied with kit.
22. Refit top section of rack using a small piece of rack to check alignment. Torque fixings to recommended 4Nm (2.95 lbf ft). Top section of rack must have at least 3 fixing screws.
23.  Refit mechanical stop.

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


27

GB

24. Thread trailing cable through trolley weight as shown.
 25. Fit trolley weight with pulley. Fit trailing cable clamp.
 - Lubricate pulley with silicone spray.
 26. Fit top end cap. Check free movement of trolley weight and cable.
 27. Secure charger and connect to fused spur/electrical connection, refer to Drawings 3109204 - Upper and 3109205 - Lower.
- Electrical Guidelines. If applicable, see Electrical drawing for details.

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GB

28. Wire charger into trailing cable plug (Ref; Wiring Diagram No. 3109208). Plug trailing cable into terminal box.
a) 0-volt
Terminate 0-volt wiring/connections as per Drawing No. 3109203.
29. Fit terminal box to underside of rail. Enclose all external wiring in flexible conduit for mechanical and weatherproof connection.
30. Fit slowing and stopping ramps to underside of rail using them to hold cables inside rail. Fit cable retaining clips as necessary.

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


**GB**

31. Fit lower section of capping strip. Fit top section, mark and cut to required length.
- a) Capping strip
32.  Adjust slowing and stopping ramps to correct positions by sliding up and down rail. Lubricate rack and spray silicone spray on trailing cable. Test unit in accordance with Test Certificate (Page 21). Refit covers and fit all labels as shown on Page 32. You are now ready to handover the lift as per Page 28.

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



NOTE :
SEE ALSO DIAGRAM 3109202
FOR FURTHER LANDING
CONTROL CONNECTION
DETAILS

NOTE :
FOR CLARITY, CHARGER CONNECTION
HAS NOT BEEN SHOWN. CHARGER
WIRES CONNECT TO TRAILING CABLE
AT CONNECTION PCB, +24V TO CON1,6,
AND 0V TO CON1,6.
SEE DRGS : 3109204
3109205

IMPORTANT NOTES :
CON3 LINK WIRE/CON8 JUMPER PLUG MUST BE REMOVED
WHEN HARD WIRED LANDING CONTROLS ARE FITTED.

PUSH ON
CONNECTIONS



NOTE :
THE 'NO CHARGE' ALARM IS SILENCED
BY THE BATTERY ISOLATOR SWITCH,
NOT THE CHAIR ARM ON/OFF OR
KEYSWITCH.

- D** Verdrahtungs- und Fehlerdiagnosedetails - 3209003
- E** Detalles de cableado y de diagnóstico de avería - 3209003
- F** Détails du câblage et de la localisation des pannes - 3209003
- I** Schemi cablaggio e risoluzione guasti - 3209003
- Jap** 配線および故障診断の詳細 - 3209003

- D Verbindungsplan für die Ladevorrichtung – obere Ebene - 3109204
- E Diagrama de conexión del cargador - Nivel superior - 3109204
- F Schéma de raccordement du chargeur – niveau supérieur - 3109204
- I Diagramma collegamento caricatore – Livello superiore - 3109204
- Jap 充電器接続図- 上部レベル-3109204


- D Verbindungsplan für die Ladevorrichtung – untere Ebene - 3109205
- E Diagrama de conexión del cargador - Nivel inferior - 3109205
- F Schéma de raccordement du chargeur – niveau inférieur - 3109205
- I Diagramma collegamento caricatore – Livello inferiore - 3109205
- Jap 充電器接続図—下部レベル— 3109205

- Verbindungsplan für die Nullspannung - 3109203
- Diagrama de conexión de voltaje cero - 3109203
- Schéma de connexion zéro volt - 3109203
- Diagramma collegamento zero Volt - 3109203
- ゼロボルト接続図-3109203

GB**Instructions (all fixings stainless steel or copper)**

- 1 Cut the 0V wire from connection PCB (1) in terminal box (2) to approximately 150 mm. Fit blue ring crimp (8) and connect wire to 0V stud (3) on underside of terminal box lid.
- 2 Drill Ø4.2 mm hole on wall side of rail (10) as shown and tap M5. Connect 0V link (6) to 0V stud (5) on upper side of terminal box lid (4). Use M5 screw (9) and externally toothed washer (7) to connect the other end of 0V link (6) to rail (10) using the drilled and tapped M5 hole (11) on wall side of rail.
- 3 Test as detailed in 310 installation manual.
 - When cables have been routed into the terminal box, the cable glands must be fully tightened to ensure a waterproof seal.
 - For clarity the terminal box connections are not shown. All connections are as 400, except for landing controls (see drawing 3109202).

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


DIRECTIONAL CONTROL ARM


SECURE UNUSED WIRES WITH CABLE TIES AS SHOWN

POWER CONTROL ARM

GB

- A** Switched fuse spur
- B** Battery charger unit
- C** Stannah installer to fit and connect
- D** ABS box secured to wall at a maximum distance of 3m from spur unit with cable protectors inside cavity space and secured internal 2 way connector block
- E** Internal
- F** External
- G** Round 2 core multi-strand to BS6004/6346 or international equivalent:
Max length of run from ABS box in 1,0mm² = 22m
Max length of run from ABS box in 1,5mm² = 30m
(External wiring shown routed to top of rail. Route to lower rail if charger at bottom of staircase)
- H** Rail terminal box
- J** Rail conduit terminal bracket.
20mm flexible conduit to be connected to rail conduit terminal bracket (fitted at lower rail position if charger at bottom of staircase)

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GB

- All internal and external electrical work shall meet the current IEE wiring regulations, British Standard BS7671. The work shall also conform to local authority specifications (where applicable), providing there is no breach of the current IEE regulations. The work shall conform to, as a minimum standard, any equivalent or additional international standard applicable for that country when installed outside of the UK.
- A NICEIC certificate is mandatory with all electrical work inside the UK. An equivalent international certificate is required for all electrical work outside of the UK.
 - All external electrical work should be protected to a minimum of IP65.
- Electrician to run length of specified 2 core cable from building to no further than 0.5m from nearest riser of staircase and leave an additional 8m coiled length for Stannah installer to connect to rail terminal box.
 - The first 3m (1) of the additional 8m (2) coiled length must be enclosed in 20mm flexible PVC conduit.

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Certificate of Test and Examination after Installation of 320 Model Stairlift to Ensure Compliance with CE Directives:

Machinery Directive (98/37/EEC) Low Voltage Directive SI (94/3260), EMC Directive SI (92/2372).

Name and address of user:..... Tel:

- 1. DESCRIPTION**
 a) Manufacturer: *Stannah Stairlifts Ltd*
 b) Stairlift Reference number:.....Chair batch No. Carriage batch No.
 c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)

- d) Rated Speed: *0.1 m/s*
 e) Motor tested on assembly YES / NO
 f) Control actuator type: *Constant pressure controls*
 g) Model Type: *300*
 h) Type of Motor: *DC*
 i) Motor protection (type): *Electronic*
 j) Type of drive: *Rack and Pinion*

2. ELECTRICAL TESTS

- a) Electrical supply voltage at time of test:..... V AC
 b) Battery voltage:..... V DC
 c) Stall current: Amps
 d) Tripping Time: Seconds
 e) Control voltage: *24V DC*
 f) Mains fuse rating: Amps
 g) Polarity check : L-N = L-E = N-E =
 h) 0-volt rail connected and continuity tested.... YES / NO
 i) Confirm battery polarity check: + ve - ve
 j) AC Supply frequency:Hz
 k) Phase: *Single*
 l) Motor reversal delay:Seconds
 m) Motor operation with rated load:

	Running current : Amps	Running speed : RPM or m/s
UP		
DOWN		

3. BATTERY CHECKS

- a) Battery charge voltage: V DC
 b) Confirm battery enclosure is vented: YES / NO
 c) Confirm battery is securely fitted: YES / NO
 d) Confirm batteries are fully charged: YES / NO
 e) Confirm no leak when battery is fitted: YES / NO
 f) Confirm off charge indicator is functional: YES/NO

4. MECHANICAL / OPERATION TEST SWITCH OPERATION

- a) Confirm all on/off stop switches operate correctly YES / NO
 b) Confirm all control limits operate correctly: upper YES / NO lower YES / NO
 c) Confirm final limits operate correctly: upper YES / NO lower YES / NO
 d) Confirm all sensitive edges/surfaces operate correctly YES / NO
 e) Confirm safety gear switch mechanism operates correctly YES / NO
 f) Confirm that all interlock switches operate correctly* YES / NO
 g) Tests to verify the function of the Overspeed governor and safety gear
 h) Undergo without failure a series of static load tests to check the forces imposed by the maximum working load

COMPONENT FUNCTIONING

- i) State the stopping distance of the carriage with rated load upon application of the brake: Up.....mm Down.....mm
 j) Confirm hand winding operates correctly YES / NO
 k) Confirm overspeed governor has been subjected to factory test for correct tripping speed: YES / NO
 l) If alarm is fitted, confirm mode of operation and that the system operates correctly: YES / NO / NA
 m) State distance of footrest above floor level: Upper.....mm Lower.....mm Intermediate.....mm

5. CONFIRMATION

Confirm operating instructions have been handed to the user: YES / NO

6. DECLARATION

7. We certify that onthis stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.
 Signed: Position/Qualifications:.....
 Address:..... Date:

Introduction

This work method is used for guidance when commissioning a stairlift and completing the Certificate of Test and Examination after the stairlift is installed.

Section 1

Description

- b) Enter details:
 • The contract or order reference number, e.g. 99390
 • The chair serial no.
 • The carriage serial no.
 e) Check that the motor tested label is present on the motor.

Section 2

Electrical Tests

- a) Measure AC supply across the live and neutral terminals at the charger supply input.
 b) Measure battery voltage.
 c) Measure the current draw from the motor during stalling.
 d) Electronic trip time (Maximum 3 seconds).
 f) Fuse rating (Max 2 Amps).
 g) Check polarity at spur or the incoming mains supply.
 h) Carry out a continuity test to ensure the 0 volt wire is mechanically and electrically connected to the rail.
 i) Check polarity of battery connections
 j) Enter mains AC supply frequency (Hz) for the appropriate country.
 k) Motor reversal delay (Min 3 seconds).
 m) Measure the current draw from the motor during operation in both directions.

Section 3

Battery Tests

- a) Measure battery charge voltage.
 b) Delete as necessary.
 c) Delete as necessary.
 d) Test battery charge condition using the four volt drop principle.
 e) Visually check for any leaks or corrosion.

Section 4

Mechanical / Operation Test

Switch Operation

- a) Delete as necessary.
 b) Delete as necessary.
 c) Delete as necessary.
 d) Delete as necessary.
 e) Confirm safety gear mechanism has been subjected to factory test (Test label present) and will function correctly when tripped.
 f) All electrical safety devices function correctly.
 g) Confirm stairlift is capable of supporting maximum rated load + 25%.
Component functioning
 h) Record the stopping distance (must not exceed 20mm).
 i) Check the handwinding operation.
 j) Confirm the tested label is present on the overspeed governor assembly.
 k) State option fitted.
 l) Record footrest heights.

Section 5

Confirmation

The stairlift user must be fully instructed on the use of the stairlift including the emergency handwinding procedures. All operating instructions should be handed over to the user on completion.

Section 6

Declaration

This section must be completed in full using the distributor's address and signed.

On completion of the final test and examination and completion of the test certificate affix the warning labels to the carriage assembly referring to the installation manual (See page 32). This is important to ensure compliance with CE directives.

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Certificate of Test and Examination after Installation of 320 Model Stairlift to Ensure Compliance with CE Directives:

Machinery Directive (98/37/EEC) Low Voltage Directive SI (94/3260), EMC Directive SI (92/2372).

Name and address of user: Tel:

1. DESCRIPTION

- a) Manufacturer: *Stannah Stairlifts Ltd*
- b) Stairlift Reference number: Chair batch No. Carriage batch No.
- c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)
- d) Rated Speed: *0.1 m/s*
- e) Motor tested on assembly YES / NO
- f) Control actuator type: *Constant pressure controls*
- g) Model Type: 300
- h) Type of Motor: *DC*
- i) Motor protection (type): *Electronic*
- j) Type of drive: *Rack and Pinion*

2. ELECTRICAL TESTS

- a) Electrical supply voltage at time of test: V AC
- b) Battery voltage: V DC
- c) Stall current: Amps
- d) Tripping Time: Seconds
- e) Control voltage: *24V DC*
- f) Mains fuse rating: Amps
- g) Polarity check : L-N = L-E = N-E =
- h) 0-volt rail connected and continuity tested.... YES / NO
- i) Confirm battery polarity check: + ve - ve
- j) AC Supply frequency: Hz
- k) Phase: *Single*
- l) Motor reversal delay: Seconds
- m) Motor operation with rated load:

	Running current : Amps	Running speed : RPM or m/s
UP		
DOWN		

3. BATTERY CHECKS

- a) Battery charge voltage: V DC
- b) Confirm battery enclosure is vented: YES / NO
- c) Confirm battery is securely fitted: YES / NO
- d) Confirm batteries are fully charged: YES / NO
- e) Confirm no leak when battery is fitted: YES / NO
- f) Confirm off charge indicator is functional: YES/NO

4. MECHANICAL / OPERATION TEST SWITCH OPERATION

- a) Confirm all on/off stop switches operate correctly YES / NO
- b) Confirm all control limits operate correctly: upper YES / NO lower YES / NO
- c) Confirm final limits operate correctly: upper YES / NO lower YES / NO
- d) Confirm all sensitive edges/surfaces operate correctly YES / NO
- e) Confirm safety gear switch mechanism operates correctly YES / NO
- f) Confirm that all interlock switches operate correctly* YES / NO
- g) Tests to verify the function of the Overspeed governor and safety gear
- h) Undergo without failure a series of static load tests to check the forces imposed by the maximum working load

COMPONENT FUNCTIONING

- i) State the stopping distance of the carriage with rated load upon application of the brake: Up.....mm Down.....mm
- j) Confirm hand winding operates correctly YES / NO
- k) Confirm overspeed governor has been subjected to factory test for correct tripping speed: YES / NO
- l) If alarm is fitted, confirm mode of operation and that the system operates correctly: YES / NO / NA
- m) State distance of footrest above floor level: Upper.....mm Lower.....mm Intermediate.....mm

5. CONFIRMATION

Confirm operating instructions have been handed to the user: YES / NO

6. DECLARATION

We certify that on this stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.

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Certificate of Test and Examination after Installation of 320 Model Stairlift to Ensure Compliance with CE Directives: Machinery Directive (98/37/EEC) Low Voltage Directive SI (94/3260), EMC Directive SI (92/2372).

Name and address of user: Tel:

- 1. DESCRIPTION**
 a) Manufacturer: *Stannah Stairlifts Ltd*
 b) Stairlift Reference number: Chair batch No. Carriage batch No.
 c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)

- g) Model Type: 300
 h) Type of Motor: DC
 i) Motor protection (type): Electronic
 j) Type of drive: Rack and Pinion

- 2. ELECTRICAL TESTS**
 a) Electrical supply voltage at time of test: V AC
 b) Battery voltage: V DC
 c) Stall current: Amps
 d) Tripping Time: Seconds
 e) Control voltage: 24V DC
 f) Mains fuse rating: Amps
 g) Polarity check : L-N = L-E = N-E =
 h) 0-volt rail connected and continuity tested.... YES / NO

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- 3. BATTERY CHECKS**
 a) Battery charge voltage: V DC
 b) Confirm battery enclosure is vented: YES / NO
 c) Confirm battery is securely fitted: YES / NO
 d) Confirm batteries are fully charged: YES / NO
 e) Confirm no leak when battery is fitted: YES / NO
 f) Confirm off charge indicator is functional: YES/NO

- 4. MECHANICAL / OPERATION TEST SWITCH OPERATION**
 a) Confirm all on/off stop switches operate correctly YES / NO
 b) Confirm all control limits operate correctly: upper YES / NO lower YES / NO
 c) Confirm final limits operate correctly: upper YES / NO lower YES / NO
 d) Confirm all sensitive edges/surfaces operate correctly YES / NO
 e) Confirm safety gear switch mechanism operates correctly YES / NO
 f) Confirm that all interlock switches operate correctly* YES / NO
 g) Tests to verify the function of the Overspeed governor and safety gear
 h) Undergo without failure a series of static load tests to check the forces imposed by the maximum working load

- COMPONENT FUNCTIONING**
 i) State the stopping distance of the carriage with rated load upon application of the brake: Up.....mm Down.....mm
 j) Confirm hand winding operates correctly YES / NO
 k) Confirm overspeed governor has been subjected to factory test for correct tripping speed: YES / NO
 l) If alarm is fitted, confirm mode of operation and that the system operates correctly: YES / NO / NA
 m) State distance of footrest above floor level: Upper.....mm Lower.....mm Intermediate.....mm

5. CONFIRMATION
 Confirm operating instructions have been handed to the user: YES / NO

6. DECLARATION
 7. We certify that onthis stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.
 Signed: Position/Qualifications:.....
 Address: Date:

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Name and address of user: Tel:

1. DESCRIPTION

- a) Manufacturer: *Stannah Stairlifts Ltd*
- b) Stairlift Reference number: Chair batch No. Carriage batch No.
- c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)
- d) Rated Speed: *0.1 m/s*
- e) Motor tested on assembly YES / NO
- f) Control actuator type: *Constant pressure controls*
- g) Model Type: *300*
- h) Type of Motor: *DC*
- i) Motor protection (type): *Electronic*
- j) Type of drive: *Rack and Pinion*

2. ELECTRICAL TESTS

- a) Electrical supply voltage at time of test: V AC
- b) Battery voltage: V DC
- c) Stall current: Amps
- d) Tripping Time: Seconds
- e) Control voltage: *24V DC*
- f) Mains fuse rating: Amps
- g) Confirm battery polarity check: + ve - ve
- h) 0-volt rail connected and continuity tested: YES / NO
- i) Confirm battery polarity check: + ve - ve
- j) AC Supply frequency: Hz
- k) Phase: *Single*
- l) Motor reversal delay: Seconds
- m) Motor operation with rated load:

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COMPONENT FUNCTIONING

- i) State the stopping distance of the carriage with rated load upon application of the brake: Up:mm Down:mm
- j) Confirm hand winding operates correctly YES / NO
- k) Confirm overspeed governor has been subjected to factory test for correct tripping speed: YES / NO
- l) If alarm is fitted, confirm mode of operation and that the system operates correctly: YES / NO / NA
- m) State distance of footrest above floor level: Upper:mm Lower:mm Intermediate:mm

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Name and address of user: Tel:

1. DESCRIPTION

- a) Manufacturer: *Stannah Stairlifts Ltd*
 b) Stairlift Reference number: Chair batch No. Carriage batch No.
 c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)

- d) Rated Speed: *0.1 m/s*
 e) Motor tested on assembly *YES / NO*
 f) Control actuator type: *Constant pressure controls*
 g) Model Type: *300*
 h) Type of Motor: *DC*
 i) Motor protection (type): *Electronic*
 j) Type of drive: *Rack and Pinion*

2. ELECTRICAL TESTS

- a) Electrical supply voltage at time of test: V AC
 b) Battery voltage: V DC
 c) Stall current: Amps
 d) Tripping Time: Seconds
 e) Control voltage: *24V DC*
 f) Mains fuse rating: Amps
 g) Polarity check : L-N = L-E = N-E =
 h) 0-volt rail connected and continuity tested.... *YES / NO*
 i) Confirm battery polarity check: + ve - ve
 j) AC Supply frequency: Hz
 k) Phase: *Single*
 l) Motor reversal delay: Seconds
 m) Motor operation with rated load:

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3. BATTERY CHECKS

- a) Battery charge voltage: V DC
 b) Confirm battery enclosure is vented: *YES / NO*
 c) Confirm battery is securely fitted: *YES / NO*
 d) Confirm batteries are fully charged: *YES / NO*
 e) Confirm no leak when battery is fitted: *YES / NO*
 f) Confirm off charge indicator is functional: *YES/NO*

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- a) Confirm all on/off stop switches operate correctly *YES / NO*
 b) Confirm all control limits operate correctly: upper *YES / NO* lower *YES / NO*
 c) Confirm final limits operate correctly: upper *YES / NO* lower *YES / NO*
 d) Confirm all sensitive edges/surfaces operate correctly *YES / NO*
 e) Confirm safety gear switch mechanism operates correctly *YES / NO*
 f) Confirm that all interlock switches operate correctly* *YES / NO*
 g) Tests to verify the function of the Overspeed governor and safety gear
 h) Undergo without failure a series of static load tests to check the forces imposed by the maximum working load

COMPONENT FUNCTIONING

- i) State the stopping distance of the carriage with rated load upon application of the brake: Up.....mm Down.....mm
 j) Confirm hand winding operates correctly *YES / NO*
 k) Confirm overspeed governor has been subjected to factory test for correct tripping speed: *YES / NO*
 l) If alarm is fitted, confirm mode of operation and that the system operates correctly: *YES / NO / NA*
 m) State distance of footrest above floor level: Upper.....mm Lower.....mm Intermediate.....mm

5. CONFIRMATION

Confirm operating instructions have been handed to the user: *YES / NO*

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GB

- Push the OSG tool into slot in the side of the carriage.
- Push bar fully down to activate OSG lever.
- Run stairlift down while holding lever down until the OSG switch cuts out stairlift. This should happen without the safety gear touching and marking the rail.
- This test checks the integrity of all components within the mechanism as well as the safety switch.
- To reset the stairlift, remove the tool and hand-wind the stairlift up until power is restored.
- Then drive the stairlift up about 500mm to completely reset the OSG.

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GB

- Demonstration.
- Explanation.
- User safety.
- User booklet.
- Paperwork completion.

1. Ask user to stand clear of the stairlift making sure that they have a clear view of the lift.
2. Demonstrate the operation of seat and footrest; advise the user to stand clear of the footrest to avoid the footrest coming into contact with the users legs.
3. Explain the L.E.D functions on the chair clearly.
4. Explain the operations and positions of all the on/off switches. Explain that the rocker switch should always be turned off before getting in/out of the chair.
5. Explain the operation and position of the directional controls on the chair.
6. Explain that if constant pressure is released from directional controls the lift will stop, and there is approximately a three second delay before the stairlift will restart.
7. Explain the swivel function (manual or powered, if fitted) and also that the lift will stop automatically at the top and bottom of the rail.
8. Ask the user to sit comfortably on the seat and explain the use of the seatbelt (the seatbelt should always be used).
9. Ask the user to operate the stairlift to the top of the stairs. Go up first so as not to block your own access.
10. Ensure that a swivel operation is carried out, and the customer alights from the lift, this will ensure that your footrest height is correct. Ensure user has turned off chair arm rocker switch.
11. Ask the user to complete a downward journey.
12. If the stairlift has a hinge then this option will have to be fully explained at this point, explain all controls, interlocks, safety features, and general operation.
13. Once user is clear on the operation of the lift show the operation and position of all safety edges.
14. Ensure that all of your instructions have been clearly understood.
15. Invite the user to sit down comfortably while you explain the user manual. Ensure that the Service/Distributor details are pointed out.
16. Give the user the top copy of the test sheet, and ask him/her to read and sign the handover sheet. Explain that the handover sheet is also a receipt for any payment collected.

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GB

Installation tools

1. 1 off 5.5 mm nutrunner.
2. 1 off 8 mm nutrunner.
3. 1 off terminal screwdriver.
4. 1 off medium pozidrive screwdriver.
5. 1 off medium flat blade screwdriver.
6. 1 off large flat blade screwdriver.
7. 1 off 3 mm T bar or Allen key.
8. 1 off 5 mm T bar or Allen key.
9. 1 off 13 mm combination spanner.
10. 1 off 10 mm combination spanner.
11. 1 off torque wrench 3/8 drive range (12 - 68 Nm) (8.8 - 50 lbf ft)
12. 1 off torque wrench 1/4 drive range (2.5 - 11.5 Nm) (1.8 - 8.5 lbf ft)
13. 1 off socket adapter ¼ - 3/8 drive.
14. 1 off ratchet wrench 3/8 drive.
15. 1 off small ratchet wrench 1/4 drive.
16. 1 off 10" extension bar.
17. 7mm Socket Bit.
18. 1 off 8 mm long series socket bit.
19. 1 off 10 mm socket bit.
20. 1 off hammer.
21. 1 pair long nose pliers.
22. 1 pair medium outside circlip pliers/universal set.
23. 1 pair side cutters.
24. 1 pair wire strippers.
25. 1 off Stanley knife (Utility Knife).
26. 1 off tap wrench.
27. 1 off M3 tap.
28. 1 off M5 tap.
29. 1 off hexagon bit set range 2.5 - 10 mm (Long Reach).
30. 1 off torx bit T25 3/8 or ¼ drive.
31. 1 off 5 mtr tape measure.
32. 1 off small spirit level.
33. 1 off 2ft approx (61cm) level with angle finder.
34. 1 off 13 mm percussion drill.
35. 1 off mains extension lead.

Special tools	Part no
1. 1 off cable stripper	620081
2. 1 off carriage link plug	4009301
3. 1 off service tool (hand controls)	3009304
4. 1 off M8x30 grub screw	501793
5. 1 off long series magnetic bit	2992266
6. 1 off digital volt meter (DVM) conforming to current IEE / BS specifications . DVM must also be capable of measuring DC current (up to 40 Amps).	N/A
7. 1 off insulation/continuity tester conforming to current IEE / BS specifications	N/A
8. 8mm Reamer	610210
9. Wago Release Pliers	610210
10. 1 off OSG activation rod	2992339

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

**GB****Important**

Ensure correct labels are fitted on completion of the commissioning phase and prior to handover to the customer.

- Determine the appropriate load label (1) according to the angle of the rail and fit label:
 - Max Load for $\leq 45^\circ$ is 135 Kg (297.6 lbs) - 10°C
 - Max Load for $\leq 52^\circ$ is 100 Kg (220 lbs) - 5°C
- Load label is always fitted to downside of carriage.
- Choose handwinding label (2) for left hand (upside) or right hand (downside) installation.
- The cover warning label (3) must be fitted to the arm as shown.
- Fit branding label (4) in recess, as shown.

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GB

33. To calculate the correct rail length (A):
- Measure the distance from point '1' to the top nosing '2' (B).
 - Add figure from table according to stair angle (C) and appropriate footrest height (D) above upper floor level. See Page 35.
34. Ensure rack is flush with edge of upper mechanical stop (3).
- Ensure rack is flush with edge of lower mechanical stop (4).
 - Drill and tap holes for upper (3) and lower (4) stops.
4.2 mm twist drill.
M5 tap.

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



GB For Distributors with unracked rails

- From the top section of rail, slide standard lengths of rack section (980mm) into the rail channel.
- The top section of rail requires a shorter length of rack, which is cut from a standard length. This last section of rack must have a minimum of three fixings per single length of rack.
- If the remaining section of rack has less than three fixing points reduce the length of the previous section to enable the final section to have a minimum of three fixing points.
- Tighten the rack fixings to a torque of 4 Nm (2.95 lbf ft).
- Cut the last section of rack which has the minimum three fixings required.
- Fit last section of rack to the rail and tighten the rack fixings to a torque of 4Nm (2.95 lbf ft).
- Remove last section of rack to fit the carriage (see page 14).
- **8 fixings as shown for a full length of rack section which uses 2 x 4 hole rack fixing strips.**
 - A 3 fixings = 1 x 3 hole rack fixing strip.
 - B 4 fixings = 1 x 4 hole rack fixing strip.
 - C 5 fixings = 1 x 2 and 1 x 3 hole rack fixing strips.
 - D 6 fixings = 2 x 3 hole rack fixing strip.
 - E 7 fixings = 1 x 3 and 1 x 4 hole rack fixing strips.
- Rack fixing strip lengths for varying lengths of racking. The minimum length of rack fixing strip is 2-hole.

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11

GB

9. Remove carriage safety pads.
 - a) The fixing on both sides of safety pad must be unclipped before lifting safety pad off carriage.
10. Fit trailing cable to carriage by threading connector through slot on underside of carriage, ensuring enough cable is passed through to reach ECU and secure with cable clip.
 - b) Route cable as shown for left hand lift.
 - c) Route cable as shown for right hand lift.

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GB

28. Wire charger into trailing cable plug (Ref; Wiring Diagram No. 3109028). Plug trailing cable into terminal box.
a) 0-volt
Terminate 0-volt wiring/connections as per Drawing No. 3109203.
29. Fit terminal box to underside of rail. Enclose all external wiring in flexible conduit for mechanical and weatherproof connection.
30. Fit slowing and stopping ramps to underside of rail using them to hold cables inside rail. Fit cable retaining clips as necessary.

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6. DECLARATION

7. We certify that onthis stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.

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Mechanical / Operation Test

Switch Operation

- a) Delete as necessary.
 b) Delete as necessary.
 c) Delete as necessary.
 d) Delete as necessary.
 e) Confirm safety gear mechanism has been subjected to factory test (Test label present) and will function correctly when tripped.
 f) All electrical safety devices function correctly.
 g) Confirm stairlift is capable of supporting maximum rated load + 25%.

Component functioning

- h) Record the stopping distance (must not exceed 20mm).
 i) Check the handwinding operation.
 j) Confirm the tested label is present on the overspeed governor assembly.
 k) State option fitted.
 l) Record footrest heights.

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Name and address of user: Tel:

1. DESCRIPTION

- a) Manufacturer: *Stannah Stairlifts Ltd*
- b) Stairlift Reference number: Chair batch No. Carriage batch No.
- c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)
- d) Rated Speed: *0.1 m/s*
- e) Motor tested on assembly YES / NO
- f) Control actuator type: *Constant pressure controls*
- g) Model Type: 300
- h) Type of Motor: DC
- i) Motor protection (type): *Electronic*
- j) Type of drive: *Rack and Pinion*

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6. DECLARATION

7. We certify that on this stairlift was thoroughly examined and found to be free from obvious defects and that the foregoing is a correct report of the results.

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Certificate of Test and Examination after Installation of 320 Model Stairlift to Ensure

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Name and address of user: Tel:

1. DESCRIPTION

- a) Manufacturer: *Stannah Stairlifts Ltd*
- b) Stairlift Reference number: Chair batch No. Carriage batch No.
- c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)
- d) Rated Speed: *0.1 m/s*
- e) Motor tested on assembly: YES / NO
- f) Control actuator type: *Constant pressure controls*
- g) Model Type: *300*
- h) Type of Motor: *DC*
- i) Motor protection (type): *Electronic*
- j) Type of drive: *Rack and Pinion*

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Name and address of user: Tel:

1. DESCRIPTION

- a) Manufacturer: *Stannah Stairlifts Ltd*
- b) Stairlift Reference number: Chair batch No. Carriage batch No.
- c) Safe working load: ≤ 42 degrees = 160kg (353lbs) ≤ 49 degrees = 135 kg (297lbs) ≤ 52 degrees = 120kg (266lbs)

- d) Rated Speed: *0.1 m/s*
- e) Motor tested on assembly *YES / NO*
- f) Control actuator type: *Constant pressure controls*
- g) Model Type: *300*
- h) Type of Motor: *DC*
- i) Motor protection (type): *Electronic*
- j) Type of drive: *Rack and Pinion*

2. ELECTRICAL TESTS

- a) Electrical supply voltage at time of test: V AC
- b) Battery voltage: V DC
- c) Stall current: Amps
- d) Tripping Time: Seconds
- e) Control voltage: *24V DC*
- f) Mains fuse rating: Amps
- g) Polarity check : L-N = L-E = N-E =
- h) 0-volt rail connected and continuity tested.... *YES / NO*
- i) Confirm battery polarity check: + ve - ve
- j) AC Supply frequency: Hz
- k) Phase: *Single*
- l) Motor reversal delay: Seconds
- m) Motor operation with rated load:

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3. BATTERY CHECKS

- a) Battery charge voltage: V DC
- b) Confirm battery enclosure is vented: *YES / NO*
- c) Confirm battery is securely fitted: *YES / NO*
- d) Confirm batteries are fully charged: *YES / NO*
- e) Confirm no leak when battery is fitted: *YES / NO*
- f) Confirm off charge indicator is functional: *YES/NO*

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- a) Confirm all on/off stop switches operate correctly *YES / NO*
- b) Confirm all control limits operate correctly: upper *YES / NO* lower *YES / NO*
- c) Confirm final limits operate correctly: upper *YES / NO* lower *YES / NO*
- d) Confirm all sensitive edges/surfaces operate correctly *YES / NO*
- e) Confirm safety gear switch mechanism operates correctly *YES / NO*
- f) Confirm that all interlock switches operate correctly* *YES / NO*
- g) Tests to verify the function of the Overspeed governor and safety gear
- h) Undergo without failure a series of static load tests to check the forces imposed by the maximum working load

COMPONENT FUNCTIONING

- i) State the stopping distance of the carriage with rated load upon application of the brake: Up.....mm Down.....mm
- j) Confirm hand winding operates correctly *YES / NO*
- k) Confirm overspeed governor has been subjected to factory test for correct tripping speed: *YES / NO*
- l) If alarm is fitted, confirm mode of operation and that the system operates correctly: *YES / NO / NA*
- m) State distance of footrest above floor level: Upper.....mm Lower.....mm Intermediate.....mm

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