


## Wireless, intelligent bed exit warning system



The patient security system is invisibly integrated into the mattress and triggers alerts wirelessly.

### More protection for patients and residents

- 360° protection
- Invisible
- Non evasive security system

### Higher care quality and security



- Short response time
- Alerts before actual bed exit
- Easy set-up
- No tripping hazard (wireless and integrated system)


### Economic benefit

- Compatible with all common paging systems
- No additional cleansing
- No derogation of mattress or bed functionality


### How KOGNIMAT works:

The alert system is invisibly integrated into the mattress. The mattress can be placed on a patient bed regardless to its orientation ①. As soon as an unwanted bed exit is detected by the mattress ② a signal is transmitted to the receiver wirelessly ③, and the paging system sets off an alert ④.







Alert when the patient or resident **is sitting on the edge of the bed (360° around the entire mattress).**




Alert when the patient or resident **is trying to climb over the bed rail.**




Alert when the patient or resident **is caught in the bed rail**






Alert when the patient or resident **is rolling out of bed.**



Alert when the patient or resident **climbs over the bed end.**

 Developed and produced in Switzerland

## KOGNIMAT is available for these mattresses:

-  Hospisoft
-  Evoba
-  Medical Basic

## Options




**Receiver LIGHT**  
Use: bedside lamp



**Receiver FLEX**  
Use: station room  
bedside lamp  
pager  
phone



**BELL**  
Use: Bell, light, TV, etc.



**DOOR- and WINDOWSENSOR**  
Use: securing doors  
and windows